

Characterization of Nuclear Materials

Methods, Instrumentation, Applications



Introduction

The world demand for electricity is projected to double by mid-century. Because of the ever-increasing price of fossil fuels and the associated environmental concerns such as carbon dioxide emissions, nuclear power generation is regaining popularity. The production of nuclear fuel is clean and relatively inexpensive compared to fossil fuels. It is in fact the only source of clean, sustainable and affordable energy which can meet current and near-term demands for electricity

generation if sufficient numbers of reactors of current and advanced design can be brought on line in a timely manner. In addition, the nuclear reactors of current design are safe and an accident such as that which occurred at Chernobyl is virtually impossible. The resurgence of nuclear power is so strong that it is being called the nuclear renaissance.


This renaissance has led to the formation of a multi-national



cooperation called the Global Nuclear Energy Partnership (GNEP). Two of the stated aims of GNEP are to:


- develop a new generation of nuclear power plants – the so-called Generation IV systems (GEN IV) in which six different reactor types are under consideration (both thermal and fast reactors), and
- reduce waste by recycling used nuclear fuel using new technologies.

The history of commercial nuclear power generation reaches back to the mid-60s. The Generation III+ reactors now under development and the GEN IV reactors under consideration are of unique design with respect to safety, performance and economics. The GEN IV reactors include the Very High Temperature Reactor (VHTR), the Sodium-cooled Fast Reactor (SFR) and perhaps the most unique, the Molten-Salt Reactor (MSR). Many of these reactors will contain unique fuel systems and reactor materials.




Time Line for Reactor Development

Introduction




source: Idaho National Laboratory



Composition of Used Fuel

The major concern with nuclear power generation today is the safe disposal of used fuels. Because of storage-related problems, recycling of used fuel is of paramount importance. Fuel recycling is not a new concept. In fact, used fuel from light-water reactors (LWR) has been reprocessed into $(\text{U,Pu})\text{O}_2$ – so-called MOX fuel – for some time by various countries. This is economically feasible because spent fuel contains not only fission products and minor actinides such as Americium, Neptunium and Curium, but a large percentage of fissile Uranium and Plutonium isotopes as well.



Closed fuel cycle with front end processing of UO_2 and back end recycling/target fabrication

It has been estimated that the effective capacity of geological repositories can be increased greatly if the long-lived minor actinides such as those mentioned above plus Pu (transuranics) are separated and fabricated into the

reprocessed fuel or targets for transmutation/consumption. Several concepts have been proposed to realize this, some of which utilize fast reactors and others thermal reactors.


Introduction

In addition, studies have shown that significant reductions in repository heat and radiotoxicity loads can be realized by placing used fuel in interim storage for a few years to allow short-lived fission products such as Cs-137 and Sr-90 to partially decay prior to separation.

Interim storage also reduces the problems associated with reprocessing fuel containing Cm, but also increases the content of high-vapor-pressure

Am-241 due to the β -decay of Pu-241.

Of course, a prerequisite for the successful design of any new reactor or fuel system as well as modernization of the existing reactor fleet is the accurate knowledge of the thermophysical properties for the materials of interest. This will necessitate the measurement of the thermophysical properties for fresh, reprocessed and used fuels as well as irradiated and unirradiated




The calculated increase in repository capacity as a function of reduced Cs, Sr, Pu, Am and Cm content with separation after 25 years
 "Separations and Transmutation Criteria...", Wigeland, et al., Nuc. Tech., 2006

reactor component materials. Property measurements on fission products and/or their surrogates, glasses, containment components and geological materials associated with long-term isolation in repositories are also of paramount importance. The properties of interest include but are not limited to the thermal conductivity, thermal diffusivity, specific heat, transformation energetics, thermal expansion, bulk density, solidus/liquidus temperatures and O/M ratio.

Clearly, measurement of these properties on the materials mentioned above will necessarily be carried out in glovebox and hot cell environments as well as in cold facilities.

Over the past few years NETZSCH-Gerätebau GmbH has become the leading supplier of thermal analysis and thermophysical properties instrumentation to the nuclear industry. This is because our instruments are reliable, robust,

accurate and easy to use. All of these attributes are necessary for instruments operating in harsh environments such as gloveboxes and hot cells. Further, the modular design of our instruments makes them ideally suited for incorporation into these environments. Finally, the costs of operating in hot cells and gloveboxes are extremely high, making downtime critical. Our exemplary global service is therefore one more reason for our success.




Transmutation Benefits of Older Fuel
 "Closed Nuclear Fuel Cycle...", Collins, et al., Atalante 2008 Int. Conf., 2008

Thermophysical Properties

Introduction

Thermophysical properties can be divided into two categories – transport and thermodynamic. Transport properties include, but are not limited to, thermal conductivity, electrical resistivity and thermal diffusivity. (Actually, thermal diffusivity is a hybrid transport/thermodynamic property.)

Thermodynamic properties include specific heat, transition energetics and thermal expansion (bulk density).



Classification of Some Thermophysical Properties

Thermal Conductivity

The thermal conductivity is perhaps the single most important thermophysical property and is paramount to the design of any system operating at elevated or sub-ambient temperatures. It consists of a lattice and/or electronic component, depending on the material (other components are also possible). It is well known in the nuclear industry that the thermal conductivity controls:

- temperature gradients in fuel
- efficiency of cladding and heat exchangers

- ability of geological repositories and container material to dissipate heat
- heat transfer in multi-layer fuel systems, e.g. TRISO.

The list is almost endless. The thermal conductivity is greatly affected by corrosion, hydriding, fouling, O/M ratio, fission product carry-over, irradiation damage, composition, porosity, etc. The thermal conductivity/thermal diffusivity of almost all nuclear materials is most efficiently measured by the laser flash technique (LFA). LFAs can readily be incorporated into gloveboxes and hot cells with the appropriate modifications.

Specific Heat and Transition Energetics

The capacity of a material to store energy is governed, in part, by the specific heat (sensible heat). It is made up of lattice, electronic and defect components, depending on the material. This property is required for design of any transient heat transfer process. It is also used to quantify surface oxidation/reduction and O/M (defects) of fuels during processing. In some cases, the specific heat can be used as an indicator of the extent of damage in post irradiation examination

(PIE), e.g. stored energy. It is also required to calculate the thermal conductivity from thermal diffusivity data.

Transition energetics (latent heat) are required to characterize solid-solid transitions, melting/solidification and decomposition. Both specific heat and transition energetics are most accurately and efficiently measured by differential scanning calorimetry (DSC). The specific heat can also be measured by the laser flash technique, albeit with reduced accuracy and only with a reduced number of data points. (With DSC, generation of a quasi-continuous set of temperature-dependent specific heat data is standard.) With the required expertise, DSCs can readily be adapted for hot work.

Solidus and Liquidus Temperatures

Solidus and liquidus temperature as well as melting temperature data are necessary to establish safe reactor operating conditions and to model accident scenarios such as loss of coolant. These temperatures are greatly affected by impurities, radiation damage, O/M ratios, burnup and, of course, composition. Surprisingly, solidus/liquidus temperatures are notoriously difficult to measure accurately. DSC is the technique most often employed for these

measurements, but care must be taken to avoid undercooling during solidification (especially critical for metal alloys). Sample time constants and temperature ramp rates must be carefully considered. Solidus/liquidus temperatures of most metal alloys can also be measured by the laser flash technique (by utilization of thermal conductivity/thermal diffusivity data) and dilatometry can be used for conductors and insulators alike. For materials that melt at ultra-high temperatures, thermal arrest is sometimes employed using optical pyrometers for the temperature measurement. All things considered, DSC is the most versatile and accurate method.

Thermal Expansion

Thermal expansion can be made up of lattice, electronic, magnetic and vacancy/interstitial components, depending upon the material and temperature. Thermal expansion data are key to both reactor and fuel design. For example, it is necessary for quantification of:

- fuel swelling during irradiation
- fuel/coating compatibility (e.g. UO_2 /graphite/SiC or ZrC)
- abrasion and corrosion coating/substrate compatibility

- densification during sintering
- thermal expansion coefficients
- volume change during melting/solidification
- bulk density

and, as stated earlier, the data can be used to determine solidus and liquidus temperatures. By far the most versatile, accurate and economical technique to measure thermal expansion is pushrod dilatometry. Dilatometers are well suited for glovebox/hot cell work.

Mass Change and Evolved Gases


Temperature-dependent mass change coupled with evolved gas analysis provides valuable information to help quantify: O/M ratio, out-gassing during fuel processing, corrosion, reduction, volatile fission products/actinides during vitrification, impurities remaining from the separation process, etc. A thermogravimetric analyzer (TGA) or simultaneous TGA-DSC (STA) instrument coupled to a quadrupole mass spectrometer (QMS), either directly or by a heated transfer line, or a TGA or STA coupled to an FTIR via a heated transfer line are widely employed for these types of analysis. As with the other techniques previously discussed, these instruments can easily be modified for hot work.

Instrumentation

Flash Technique

The flash method is the fastest and most accurate way of measuring the thermal diffusivity (\Rightarrow thermal conductivity). It has been estimated that 80-85% of all thermal diffusivity/thermal conductivity measurements are carried out using this technique. This is because of:

- easy sample preparation: coin-sized disk or square/rectangle (can be even smaller)
- flexibility: measurements on solids, liquids, pastes, powder and laminate samples
- small sample size: ideal for fuels/irradiated materials
- fast measurement times: due to small sample size
- high accuracy: usually better than 5%
- wide thermal diffusivity/thermal conductivity range: approx. 0.01 to 1000 W/m·K
- wide temperature range: -125°C to 2800°C



Comparison of the flash method with other techniques utilized in NETZSCH instruments as well as approximate room temperature thermal conductivity ranges of several types of materials.



As shown below, the front surface of a sample is heated by a short energy pulse (xenon or laser) and the time-dependent temperature rise on the rear surface is measured by an IR detector. From a plot of T/T_{\max} (or V/V_{\max}) vs half-time, the thermal diffusivity of a perfectly insulated sample at one half-time is given by:

$$a = 0.1388 l^2/t_{1/2}$$

where

- a = thermal diffusivity (mm^2/s)
- l = sample thickness (mm)
- $t_{1/2}$ = half-time (s)


There are several sophisticated models to correct for heat loss (non-adiabatic conditions) and the finite pulse width of the energy source.

Using the measured thermal diffusivity values, the thermal conductivity is calculated by:

$$\lambda(T) = a(T) \rho(T) c_p(T)$$


where

- λ = thermal conductivity (W/m-K)
- ρ = bulk density (g/m^3)
- c_p = specific heat (J/g-K)



Experimental Data

The bulk density can be obtained from thermal expansion data and the specific heat by either the LFA or a DSC.



Flash Technique

Instrumentation


NETZSCH LFAs

LFA 457 *MicroFlash*[®]

The 457 embodies the latest technology for modern laser flash systems. This bench-top instrument has a temperature range of -125°C to 1100°C using two user-exchangeable furnaces. The vacuum-tight design enables tests under well-defined atmospheres. Sample sizes of up to 25.4 mm Ø can be accommodated and the integrated sample changer allows simultaneous measurement of up to 3 samples. The innovative infrared sensor technology employed in the system enables measurement of the sample temperature increase, even at temperatures of -125°C.

The vertical arrangement of the sample holder, furnace and detector simplifies sample placement and, at the same time, guarantees an optimum signal-to-noise ratio for the detector signal. The small footprint of this instrument makes it ideal for modification for glovebox and hot cell applications.

LFA 457 modified for glovebox operation



LFA 447 *NanoFlash*[®]

The 447 is a compact xenon flash-based system with a temperature range of RT to 300°C. The temperature rise is measured by an IR detector. The integrated automatic sample changer allows unattended analysis of up to 4 samples. The 447 accommodates samples with multiple geometries and sizes up to 25.4 mm Ø.

The optional MTX (scanning) extension allows high-resolution measurement of the spatial distribution of the thermal diffusivity/conductivity on large 5 cm × 5 cm samples.

The NETZSCH LFAs cover almost the entire range of temperatures and materials of interest to the nuclear industry. For convenience, some of the specifications are summarized in the table below.

Model	Temperature Range	Atmosphere	Energy	Detector
447	RT to 300°C	inert, oxid.	Xenon flash	In-Sb
457	-125°C to 1100°C	inert, oxid., red., vac., corr.	Laser	In-Sb / MCT
427	-120°C to 2800°C	inert, oxid., red., vac., corr.	Laser	In-Sb / MCT

LFA 427

The 427 is the most powerful and versatile LFA system for research and development as well as all applications involving the characterization of standard and nuclear materials.

The LFA 427 guarantees high precision and reproducibility, short measurement times and defined atmospheres over a temperature range of -120°C to 2800°C. Special holders for liquid, fiber, paste, powder and laminate samples are available. Even fuel fragments can be tested.

The system is vacuum tight to 10^{-5} mbar. The variable laser power and pulse width functions make it easy to optimize test parameters.

LFA 427 modified for hot cell operation

Additional Information

visit www.lfa427.info



Instrumentation

Differential Scanning Calorimetry (DSC)

Heat flux DSCs function by measuring the temperature difference (voltage difference) between a sample and a reference while being heated or cooled at a constant rate. The resulting signal is proportional to the mass × specific heat (sensible heat) plus the energy associated with any phase change, etc. (latent heat). This signal is converted to actual physical quantities through calibration using melting point standards having well-defined heats of fusion or specific heat standards with well-established and reproducible temperature-dependent specific heat values. This is typically referred to as a sensitivity calibration. The signal generated by a heat flux DSC is shown here.

The temperature-dependent specific heat is obtained by running a baseline (empty crucibles), standard and sample scan. Using these data the specific heat is calculated by the so-called ratio method where the specific heat of the sample is given by:


$$C_{p,s} = C_{p,std} \frac{m_{std}}{m_s} \frac{\mu V_s - \mu V_{BL}}{\mu V_{std} - \mu V_{BL}}$$

where

$C_{p,s}$ = sample specific heat (J/g-k)
 $C_{p,std}$ = standard specific heat (J/g-k)
 m_s = sample mass (g)
 m_{std} = standard mass (g)
 μV_s = sample scan signal
 μV_{std} = standard scan signal
 μV_{BL} = baseline scan signal

The resulting specific heat curve will also include any latent heat due to phase transitions, etc. In the temperature

region of a phase transition the result is typically referred to as apparent specific heat. The resulting peak in the specific heat curve can then be integrated to give the transition energy, Δh (J/g). Finally, the onset and end of the latent heat peaks can be used to define the solidus and liquidus temperatures as long as time constants and undercooling are given proper consideration.



Heat Flux DSC Signal

The NETZSCH F Series DSCs cover most applications of interest to the nuclear industry. An overview of some of the specifications is given in this table.

Model	Temperature Range	Atmosphere	Vacuum	Sensor Types
404 F1	-150°C to 2400°C	inert, oxid., red., vac., corr., humid	10 ⁻⁴ mbar	E, K, S, B, W-Re
404 F3	-150°C to 2400°C	inert, oxid., red., vac., corr., humid	10 ⁻² mbar	E, K, S, B, W-Re

NETZSCH DSC 404 F Pegasus Series

The **F** Series DSCs (**F1** and **F3**) have a temperature range of -150°C to 2400°C and are vacuum tight down to 10⁻⁴ mbar. Eight different furnaces are available and, when equipped with a double hoist, 2 furnaces can be mounted simultaneously. There are 13 different user-exchangeable fast-fix plug-in sensors which allow optimization of sensitivity and time constants.

The **F** Series was designed specifically for high-temperature specific heat measurements. Accuracies in the range of ±2.5% can be obtained for most materials over the temperature range of -150°C to 1400°C for specific heat and better than 2.0% for latent heat over the temperature range of -150°C to 1650°C.

Evolved gases can also be identified and quantified (with optional *PulseTA*® unit) by coupling the DSC to a QMS and/or an FTIR via a heated furnace adapter head and heated transfer line. Special options for corrosive gases are also available.

Automatic sample changers are available for the **F** Series, but are not generally recommended for glovebox/hot cell operation.

These units are widely used in the nuclear industry in both cold and hot applications. The glovebox/hot cell versions are equipped with special sample loading mechanisms and sample carrier adjustment systems to facilitate operation.



DSC 404 **F1** Pegasus® with automatic sample changer

Additional Information

visit www.dsc404.info

Dilatometry (DIL)

Pushrod dilatometers measure, via an LVDT, the thermal expansion of a sample while it is being heated at a constant rate or held isothermally (with the exception of non-linear rates in rate-controlled sintering studies). The voltage signal is proportional to the expansion of the sample. This raw signal is converted to the length change by calibrating the instrument using a standard with a well-defined and reproducible thermal expansion.

The temperature-dependent length change is used to calculate the coefficient of linear expansion, α , i.e.:

$$\alpha = \frac{1}{L_0} \frac{\Delta L}{\Delta T} = \frac{1}{L_0} \frac{(L_2 - L_1)}{(T_2 - T_1)}$$

where

L = length (mm)
L₀ = initial length (mm)
T = temperature (°C)

The volumetric expansion for an isotropic material is then calculated by:

$$\frac{\Delta V}{V_0} \approx 3\alpha\Delta T$$

where

V = volume (cm³)
V₀ = initial volume (cm³)

Thus, the instantaneous temperature-dependent bulk density, ρ , for an isotropic material can be calculated using:

$$\rho = \frac{m}{V} \approx \frac{m}{V_0 (1 + 3\alpha\Delta T)}$$

where

m = instantaneous mass (g).


Additional Information

visit www.dil402.info



NETZSCH DIL 402 Series

NETZSCH manufactures 5 unique horizontal dilatometer models which operate over a temperature range of -260°C to 2800°C. The LVDT housings are thermostatted and have Invar LVDT mounting systems. This ensures that the length signal is not disturbed by the furnace or fluctuations in the room, glovebox or hot cell temperature. All of these dilatometers (except the ED) can be equipped with at least 2 user-exchangeable furnaces and some up to 4. All models except the PC are vacuum tight down to 10^{-4} mbar or better. The accuracy is better than 1% for most materials. Graphite, Al_2O_3 and fused silica sample carriers and pushrods with sample thermocouple types E, S, B and W-Re are available. Special holders for liquid metals, pastes, powders, fibers and thin films are also available.



Dilatometer Model 402 C modified for glovebox operation

Software options such as *c-DTA*[®], RCS, density and thermokinetics allow calculation of phase transitions, measurement of densification under variable heating rate conditions, calculation of bulk density from linear expansion data for isotropic samples and kinetic evaluation during phase

transitions and sintering, respectively. Using special sample holders, the solidus and liquidus temperatures can also be determined. The small footprint and horizontal design of the NETZSCH dilatometers make them ideal for incorporation into gloveboxes/hot cells.

The NETZSCH dilatometers cover almost all conceivable nuclear applications. For convenience, some of the specifications are summarized in this table.

Model	Temperature Range	Atmosphere	Δl Resolution	Type
DIL 402 PC	RT to 1600°C	inert, oxid.	8 nm	single pushrod
DIL 402 C	-180°C to 2000°C	inert, oxid., red., vac., corr., humid	0.125 / 1.25 nm	single pushrod
DIL 402 E	-180°C to 2800°C	inert, oxid., red., vac., corr., humid	1.25 nm	single pushrod
DIL 402 CD	-180°C to 2000°C	inert, oxid., red., vac., corr., humid	0.125 / 1.25 nm	dual sample/differential
DIL 402 ED	-260°C to 30°C	inert, oxid., red., vac., corr., humid	0.125 / 1.25 nm	dual sample/differential

Instrumentation

Thermogravimetric Analysis (TGA)/Evolved Gas Analysis (EGA)

TGAs measure the mass change of a sample during heating/cooling or in an isothermal phase via a sensitive microbalance. In systems in which the balance is located on top, the sample crucible is suspended in the furnace by a wire or rod. Instruments with such an arrangement are generally referred to as “hang-down” systems. If the balance is located below the furnace, the sample


crucible sits on a sensor at the terminal end of a ceramic rod. These systems are commonly referred to as “top-loading”. Some instruments make use of a horizontal design.

When TGA is coupled with DSC, simultaneous measurement of mass change and energetics (simultaneous thermal analysis - STA) with one

instrument is possible. The advantage of this is that both the mass change and energetics are measured simultaneously on one sample so that the test conditions, e.g. gas flow rate, gas partial pressure, heating rate, temperature, etc. are identical. The rate of mass change is calculated from the mass signal.

If the STA is coupled to a QMS and/or an FTIR, the mass change, energetics and evolved gases are measured simultaneously. Of course, this is ideal because the TGA and DSC tell us how much, while the evolved gases tell us what and, generally, why.

STAs can also be run in the so-called rate controlled mass change mode. In this case, the sample is not heated at a constant rate, but in such a fashion as to produce a constant rate of mass change. This is ideal for cases where high rates of mass loss can cause sample damage and/or skew the results, e.g. binder burnout during sintering.



Measured and Calculated STA Curves

Additional Information

visit www.ega-ms-ftir.info
and www.sta449.info

STA 449 **F1 Jupiter**[®]
coupled in parallel with a
QMS 403 C **Aéolos**[®] and
Bruker Tensor 27 FTIR



The NETZSCH STA/EGA systems handle almost all nuclear applications. For convenience, some system specifications are summarized below.

Model	Temperature Range	Atmosphere	Digital Resolution	MS Range
449 F1	-150°C to 2400°C	inert, oxid., red., vac., corr., humid	0.025 µg	1 - 300 amu
449 F3	-150°C to 2400°C	inert, oxid., red., vac., corr., humid	1.0 µg	1 - 300 amu
409 CD	-180°C to 2400°C	inert, oxid., red., vac., corr., humid	< 2.0 µg	up to 1024 amu
429 CD	-180°C to 2400°C	inert, oxid., red., vac., corr., humid	< 5.0 µg	up to 1024 amu

NETZSCH STA/EGA Systems

NETZSCH produces 4 different STA systems which cover a temperature range of -180°C to 2400°C. There are 9 unique furnaces available and 16 different user-exchangeable fast-fix plug-in sensors

which allow sensitivity/time constant optimization. The balances are all top-loading and purge gases flow from the bottom and exit at the top of the furnace along with the evolved gases.

This arrangement is ideal for EGA work (chimney effect). The balances have high capacities and digital resolutions (see above). The units are all vacuum tight, making them well suited for work requiring controlled atmospheres. The STAs have excellent long-term stability (low drift) because the balance chambers are thermostatted to isolate the balance from fluctuations in glovebox/hot cell temperatures. The vacuum-tight design and long-term stability make the NETZSCH STAs ideal for corrosion and O/M studies. An optional water vapor furnace and generation system is an additional advantage for corrosion studies.

Finally, the NETZSCH STAs can be coupled with MS and/or FTIR via a heated adaptor head and heated capillary. The STA 409 CD and STA 429 CD can be directly coupled to an MS via Skimmer® technology. A well-defined quantity of gas can be injected into the STA via a *PulseTA*® unit. This permits calibration of the MS or FTIR so that the gases released from the sample can not only be identified, but quantified as well.




Applications

Thermal Conductivity

Thermal conductivity of four stoichiometric oxide ceramics


Depicted here is the thermal conductivity of four stoichiometric oxide ceramics with 95% to 98% theoretical density. The materials show the classical trends for non (or low) electrical conductors. That is, exponential and 1/T decay with increasing temperature due to defect and phonon-phonon scattering, respectively. Further, decreasing thermal conductivity with increasing difference in ion mass is typical, that is:

$$\lambda_{\text{UO}_2} < \lambda_{\text{ThO}_2} < \lambda_{\text{Al}_2\text{O}_3} < \lambda_{\text{BeO}}$$




Thermal conductivity of some actinide nitrides

Increasing thermal conductivity with temperature indicates a significant electronic component and decreasing electrical resistivity. Decreasing thermal conductivity with increasing ion mass difference suggests increased phonon-electron scattering, i.e. shorter mean free path.



Thermal conductivity of three carbide ceramics


Even though a conductor, SiC follows the classical lattice thermal conductivity trend, while ZrC displays more of a dominant electronic trend. UC shows only a slight decrease in thermal conductivity with increasing temperature. This is due to the fact that the electronic and lattice component trends almost nullify each other even though the electronic component is dominant.



Thermal conductivity of Uranium Nitride

In this example, the thermal conductivity of UN has been separated into the lattice and electronic components. Here the typical lattice 1/T and increasing electronic behavior is displayed. Note the noncoincidental agreement between the lattice thermal conductivities of UN and UO₂. It must be pointed out here that

the lattice and electronic components of thermal conductivity were calculated from thermal conductivity and electrical resistivity measurements using a fitting routine. The Wiedemann-Franz-Lorenz (WFL) law, which can lead to significant errors in regions of non-elastic scattering, was not applied directly.




Applications

Specific Heat and Transition Energetics

Specific heat of graphite


The specific heat of graphite was measured over the temperature range of -100°C to 1300°C . As expected, the trend follows the Debye theory (electronic specific heat is negligible at these temperatures). The measured and published values deviate by a maximum of $\approx 5.0\%$ and are generally within $\approx 2.0\%$. Note the reproducibility of the measurements over both the low-temperature (LT) and high-temperature (HT) ranges.



Heat flow in a zirconium alloy

The heat flow during heating and cooling for a Zr-Nb alloy over the temperature range of 400°C to $\approx 1050^{\circ}\text{C}$ has been analyzed. Of interest here is the solid-solid endothermal phase transition between $\approx 815^{\circ}\text{C}$ and 990°C on heating. Upon cooling, the


exothermal phase transition is shifted to lower temperatures, i.e. between $\approx 985^{\circ}\text{C}$ and 750°C . Notice that the transition energetics are relatively equal at 47.4 J/g and 47.5 J/g on heating and cooling, respectively. This is not atypical for metal alloys.



Thermal Expansion


Thermal expansion of POCO graphite

The thermal expansion of high-purity POCO graphite was measured over the temperature range of 500°C to 2800°C. Clearly the reproducibility is excellent and, as shown, the deviation from literature values is significantly less than 1.0%.



Thermal expansion of some nitride ceramics

In comparing the thermal expansion of several nitride ceramics, it is interesting to note that the thermal expansion increases with ion mass difference. This increasing asymmetric behavior is also responsible for increased phonon and phonon-electron scattering, resulting in shorter mean free paths and therefore inverse behavior between thermal expansion and thermal conductivity in many ceramics.





Applications

Urania


Mass change of hyperstoichiometric UO_2

Plotted are the mass change and rate of mass change in inert and reducing atmospheres. Clearly the curves overlay until $\approx 325^\circ\text{C}$. The first mass loss step below 325°C is due to the evolution of moisture and impurities. Above 325°C the curves deviate because of the greater oxygen loss in the reducing atmosphere. By 700°C the urania is stoichiometric in the reducing atmosphere, while in the inert atmosphere the material remains over-stoichiometric. This has far-reaching consequences with regard to the sinterability.



Oxidation of stoichiometric UO_2 in an air atmosphere

The conversion from UO_2 to U_3O_7 to U_3O_8 is of course accompanied by mass increases of 1.14% and 1.72%. The exothermal energetics of the transitions are clearly shown by the DSC curve. It is interesting to note the excellent correlation between the mass change, rate of mass change and energetic curves. This plot clearly shows the advantage of simultaneous TG-DSC (STA) measurements using one instrument.




Waste Surrogates

Thermal diffusivity and specific heat of a waste glass

The thermal diffusivity and specific heat of a waste glass consisting of oxides ranging from alumina to neodymia are plotted here. The specific heat data show the glass transition (endothermic) in the range of 500°C followed by cold crystallization (exothermic) between 600°C and 750°C, indicating that the sample was cooled fast enough to produce a significant amorphous phase. Also, the large excess enthalpy peak superimposed on the glass transition may indicate significant aging time below the transition. The diffusivity displays the classical 1/T behavior followed by a significant drop across

the glass transition. The correlation between the two data sets is excellent.




Evolved gas measurement on surrogate waste material

Waste material can contain silicates like clays contaminated with radioactive inorganic salts (hydrates, sulphates, nitrates etc.). The processing of these waste materials is often done by a heat treatment. TG-MS measurements allow the determination, identification and quantification of the evolved gases during heating of these materials. The surrogate material for this investigation contained vermiculite and alkali and earth alkali nitrates. During linear heating at 10 K/min in a He atmosphere, water, nitrogen oxide, nitrogen dioxide and carbon dioxide were detected. The water came off between RT and about 400°C, nitrogen dioxide shows two intensity maxima at about 180°C and 497°C. The nitrogen oxide emission starts at around 300°C and


shows maximum intensity at 658°C. The amount of CO₂ is very small. Most probably there are some impurities of organic material, hydrocarbonates or

carbonates. Even during storage in an air atmosphere some carbonates could be formed due to the small grain size of the material.




Applications

Ceria



Densification of CeO₂

The shrinkage and rate of shrinkage of ceria at five different heating rates during sintering are plotted here. The samples were pressed to 51-53% theoretical density at 355 MPa. No binder or mold lubricants were employed. The onset of densification occurs at ≈800°C and is not complete by 1660°C. Typical heating rate dependence due to rate-limited metal diffusion is displayed. The curves show a two-step densification process which is somewhat atypical. Final density was 81-85% of theoretical.




O/M ratio of CeO₂ during sintering

This figure shows the O/M ratio during heating. These values were calculated from TGA data measured under multiple partial pressures of oxygen (PO₂). Sample preparation was identical to that above. During heating the O/M starts to decrease at ≈1000°C and there are clearly different reduction rates resulting from the variable PO₂ over the sample. The process is reversible. The resulting deviation from stoichiometry has a profound effect on the sinterability.


Specific heat of CeO₂ during sintering

Depicted here is the specific heat of ceria with and without reduction. Sample preparation remained unchanged. In the absence of reduction the specific heat follows the Debye theory and is in good agreement with published values. With reduction the specific heat (actually apparent specific heat) increases non-linearly with temperature. This is a result of the energy required for the increased rate of defect formation. The reduction becomes obvious at $\approx 950^\circ\text{C}$. The difference between the two curves is often referred to as the excess specific heat.



Property evolution of CeO₂ during sintering

In this figure the thermal diffusivity, specific heat, bulk density, rate of change of bulk density and O/M are superimposed. First, the correlation between the specific heat and O/M is obvious. Next, the thermal diffusivity increases between $\approx 400^\circ\text{C}$ and 850°C due to surface diffusion (no densification). The increase in thermal diffusivity above $\approx 850^\circ\text{C}$ is a result of decreasing porosity (densification), but is partially attenuated by the decreasing O/M (increased defect scattering). The densification (diffusion) rate between $\text{O/M}=1.9995$ and 1.9840 (1165°C and 1620°C) is controlled, in part, by the O/M-created defects. Likely controlling mechanisms are $V_M \rightarrow M_i$ cluster. XRD measurements show that



above $\approx 1600^\circ\text{C}$ coarsening occurs at the expense of densification, which manifests itself as a drop in the sintering rate.

Thermophysical Properties of Some Important Nuclear Materials at Room Temperature

The properties given in this table are intended to be used as estimates only. Some are quite accurate while others are less certain. Of course it is well known that small changes in stoichiometry, porosity, microstructure, etc. can have a profound effect on the thermophysical properties. As a result, NETZSCH does not guarantee the accuracy of these values.

Fuels			
Material	Property		
	λ (W/m·K)	C_p (J/g·K)	ρ (g/cm ³)
Pu	6.7	0.145	19.816
PuC	6.9	0.176	-
PuN	9.6	-	14.400
PuO ₂	6.3	0.240	11.460 th.
Th	54.0	0.118	11.700
ThO ₂	13.0	0.235	9.110 th.
U	27.6	0.116	19.070
U-14%Zr	18.5	-	-
U-91%Zr	7.4	-	-
(U _{0.8} Pu _{0.2})O ₂	8.9 @ th.	0.240	11.660 th.
UC	25.3	0.200	13.630 th.
UN	13.0	0.190	14.300 th.
UO ₂	8.7 @ th.	0.235	10.960 th.

Misc.			
Material	Property		
	λ (W/m·K)	C_p (J/g·K)	ρ (g/cm ³)
Al ₂ O ₃	36.0	0.765	3.970
MgO	60.0	0.929	3.770
Si ₃ N ₄	16.0	0.691	2.700
SiC	120.0	0.675	3.160
ZrC	26.3	0.378	6.730
ZrN	10.5 (80-90% th.)	0.386	7.090
ZrO ₂ (Y-stab.)	1.9	0.455	5.680

Structural			
Material	Property		
	λ (W/m·K)	C_p (J/g·K)	ρ (g/cm ³)
InX-750	11.7	0.439	8.510
Low Cr Steel	42.3	0.442	7.858
Nichrome	12.0	0.420	8.400
SS-304	14.9	0.477	7.900
SS-316	13.4	0.468	8.238
Zirc2	17.0	0.290	6.550
Zirc4	14.1	0.293	6.580
Zr-2.5%Nb	19.0	0.315	6.570

Additional elements

Material	Property		
	λ (W/m-K)	C_p (J/g-K)	ρ (g/cm ³)
Ag	429.0	0.235	10.500
Al	237.0	0.903	2.702
Cd	96.6	0.231	8.650
Ce	11.3	0.192	6.770
Cr	93.7	0.449	7.160
Fe	80.2	0.447	7.870
In	81.8	0.267	7.310
Mg	156.0	1.024	1.740
Mo	138.0	0.251	10.240
Nb	53.7	0.265	8.570
Ni	90.7	0.444	8.900
V	30.7	0.489	6.100
W	174.0	0.132	19.300
Zr	22.7	0.278	6.570


Moderators

Material	Property		
	λ (W/m-K)	C_p (J/g-K)	ρ (g/cm ³)
Be	200.0	1.825	1.850
BeO	270.0	1.030	3.000
D ₂ O	0.560	4.230	1.100
Graphite	100.0	0.710	1.700
H ₂ O	0.561	4.217	1.000
ZrH _{1.8}	34.0	0.410	5.620

Absorbers

Material	Property		
	λ (W/m-K)	C_p (J/g-K)	ρ (g/cm ³)
AgInCd	60.0	0.230	10.170
B	27.0	1.107	2.500
B ₄ C	92.0	0.960	2.520
BN	28.7	0.848	2.250
Hf	22.3	0.363	13.090
HfB ₂	22.6	0.396	11.200
Ta	57.5	0.140	16.600
ZrB ₂	23.0	0.230	6.090

Expertise in Service



Our Expertise – Service

All over the world, the name NETZSCH stands for comprehensive support and expert, reliable service, before and after sale. Our qualified personnel from the technical service and application departments are always available for consultation.

In special training programs tailored for you and your employees, you will learn to tap the full potential of your instrument.

To maintain and protect your investment, you will be accompanied by our experienced service team over the entire life span of your instrument.

Summary of Our Services

- Installation and commissioning
- Hotline service
- Preventive maintenance
- Calibration service
- IQ /OQ/PQ
- On-site repairs with emergency service for NETZSCH components
- Moving/exchange service
- Technical information service
- Spare parts assistance

Our Expertise – Applications Laboratories

The NETZSCH Thermal Analysis Applications Laboratories are a proficient partner for nearly any thermal analysis issue. Our involvement in your projects begins with proper sample preparation and continues through meticulous examination and interpretation of the measurement results. Our diverse methods and over 30 different state-of-the-art measuring stations will provide ready-made solutions for all your thermal needs.

Within the realm of thermal analysis and the measurement of thermophysical properties, we offer you a comprehensive line of the most diverse analysis techniques for materials characterization (solids, powders and liquids).

Measurements can be carried out on samples of the most varied of geometries and configurations. You will receive high-precision measurement results and valuable interpretations from us in the shortest possible time. This will enable you to precisely characterize new materials and components before actual deployment, minimize risks of failure, and gain decisive advantages over your competitors.

For production problems, we can work with you to analyze concerns and develop solutions. The minimal investment in our testing and services will reward you with reduced down time and reject rates, helping you optimize your processes across the board.



The NETZSCH Group is an owner-managed, internationally operating technology company headquartered in Germany.

The three Business Units – Analyzing & Testing, Grinding & Dispersing and Pumps & Systems – provide tailored solutions for highest-level needs. Over 2,300 employees at 130 sales and production centers in 23 countries across the globe guarantee that expert service is never far from our customers.

Leading Thermal Analysis

When it comes to Thermal Analysis, Adiabatic Reaction Calorimetry and the determination of Thermophysical Properties, NETZSCH leaves no stone unturned. Our 50 years of applications experience, broad state-of-the-art product line and comprehensive service offerings ensure that our solutions will not only meet your every requirement but will also exceed your every expectation.

NETZSCH-Gerätebau GmbH
Wittelsbacherstraße 42
95100 Selb
Germany
Tel.: +49 9287 881-0
Fax: +49 9287 881 505
at@netsch.com

www.netsch.com