

Leading Thermal Analysis.

NETZSCH

NEW

Proteus® Software Version 4.8.5*

Update

- now also available for Windows Vista™ Business and Windows Vista™ Enterprise

Get in touch with us!

*) for users of Proteus Versions 4.XXX

NETZSCH-Gerätebau GmbH
Wittelsbacherstraße 42
95100 Selb
Germany

Telephone: +49 9287 881 110
Fax: +49 9287 881 505
E-Mail: at@netsch.com

www.netsch.com

Features

Proteus® Software Version 4.8.5

- Determination of peak and onset temperatures as well as other temperature values, and depiction of these on a time scale
- Free selection of the scaling of the DSC integral curve as a % or as J/g
- Onset evaluations on the temperature curve (e.g. for Photo-DSC)
- SuperRes® now also for the TG 209 F3 Tarsus
- ... and other useful functions

Also note the already existing features from previous 4.8.x versions such as:

- USB support also for older TASC 414/2/3/4 instruments (adapter available)
- ASC measuring software with new, flexible features:
 - The ASC program can be fully edited during a running measurement
 - Pauses can be programmed
 - The programmed sequence can be altered while in progress; new measurements can be inserted and/or the set order can be changed
 - ...

An additional advantage of the Proteus® Software:
You can simply install the software yourself.
No costly visit by a service technician is necessary.