


## Thermoanalytical Characterization of Pharmaceuticals (engl.)

NETZSCH Jahrbuch 2000



The NETZSCH Annual 2000 **Thermoanalytical Characterization of Pharmaceuticals** (engl.) is based on decades of experience of the renowned specialists Dr. E. Marti, Solvias AG (former Novartis), Basle / Switzerland, and Prof. Dr. W.-Y. Ma, Shanghai Institute of Pharmaceutical Industry, Shanghai / P.R.C., in the fields of Pharmaceutical Research and Technology, as well as E. Kaisersberger and Dr. G. Kaiser, Netzsch "Thermal Analysis Applications", NETZSCH-Gerätebau GmbH, Selb / Germany.

Selected examples of drug substances and excipients are outlined in this workbook, which is based on thermoanalytical measurements using DSC, thermobalance, TG-FTIR, TG-DSC-MS and c-DTA.

Each substance is introduced by literature data and pharmaceutical use. Special attention having been paid to the presentation of easily understandable experimental results.

In addition, the collected experimental facts have been evaluated by different software programs. Routines have been applied for the determination of melting points, enthalpy of fusion and crystallization, dehydration, purity, kinetic data, Gibbs free energy function, melting range and molten fraction. The importance of gas analysis by FTIR and mass spectrometry has been underlined.

The results obtained with the performed investigations are brought into the context of thermodynamics, physical chemistry and solid-state or material sciences.

Below you can have a look at the [table of contents](#)

Bestellen Sie noch heute Ihr Exemplar  
per E-mail: [at@netzsch.com](mailto:at@netzsch.com)  
per Fax: (+49) 9287 / 881-144  
Preis: 49,- Euro zzgl. MWSt. und Versandkosten

# Thermoanalytical Characterization of Pharmaceuticals

## Table of contents

### Applications:

No.	Substance	Topics
1.	Cimetidine A	Melting point Glass transition
2.	Paracetamol Polymorphism	Thermal resistance
3.	Carbamazepine	Polymorphism Glass transition Gibbs free energy functions Stability of polymorphs Thermal stability
4.	Sulfathiazole	Melting, polymorphic Transformation Gibbs free energy functions Stability of polymorphic forms
5.	Polyethylene Glycol 6000	Melting point Liquid fraction
6.	Polyethylene Glycol 1500	Melting point Liquid fraction
7.	Microwax-Paraffin	Melting range Liquid fraction Multicomponent Analysis
8.	Naftopidil	Eutectic purity, solvent content
9.	Zaleplon	Eutectic purity in open and closed Sample pans, Solvent content Evaporation of solvent in closed pan
10.	Vanillin Ribavirin Sulfadimidin Saccharin Phenolphthalein	Melting point, melting range: comparison between DSC and Capillary method
11.	Phenacetin p-Aminobenzoic Acid	Eutectic phase diagram: Determination by DSC and Theoretical approach
12.	Castor Oil	Thermal, oxidative stability
13.	Calcium-Oxalate Monohydrate	Mass calibration Decomposition process by TG-FTIR
14.	a-D-Glucose	Thermal decomposition by DSC, TG-DSC-MS, Kinetic evaluation, Predictions
15.	Estradiol Hemihydrate	Dehydration Melting, c-DTA
16.	Ranitidine Hydrochloride	Thermal stability Melting point: function of the heating rate Decomposition by TG-FTIR Glass transition
17.	a-Lactose Monohydrate	Thermal characterization for Hydrate, anhydrate DSC, TG-FTIR
18.	Amorphous a-Lactose e	Glass transition and Amorphicity for the Monohydrate and anhydrat

FAX Nr.: (+49) 9287 / 881 144



## Thermoanalytical Characterization of Pharmaceuticals

Ich bestelle ..... Exemplar/e

**Thermoanalytical Characterization of Pharmaceuticals**

jedes für 49,- Euro, zzgl. MwSt. und Versandkosten.

Ich erhalte eine Rechnung mit Zustellung der Ware.

Meine Anschrift:

Firma: \_\_\_\_\_

Name: \_\_\_\_\_ Abteilung: \_\_\_\_\_

Straße: \_\_\_\_\_

PLZ/Ort: \_\_\_\_\_

Land: \_\_\_\_\_

Telefon: \_\_\_\_\_ Fax: \_\_\_\_\_

E-mail: \_\_\_\_\_

VAT-Nr.: \_\_\_\_\_

Datum: \_\_\_\_\_ Unterschrift: \_\_\_\_\_